RECOMMENDATION FOR SANITIZATION OF COCHLEAR™ PRODUCTS IN CLASSROOMS



MINI MIC 2+

These recommendations apply to the following Cochlear™ products when shared between users:

Mini Mic 2+

A surface disinfectant such as a solution with 60-70% isopropanol should be used to disinfect the Mini Mic 2+.

When a student passes their Mini Mic 2+ to their teacher, the mic should be sanitized.

- The teacher should clean their hands with hand sanitizer prior to starting the sanitization process.
- Spray a piece of paper towel with the isopropanol disinfectant.
- Wipe the device off with the disposable paper towel.
- Never directly spray disinfectant directly onto the device.
- Make sure openings, in particular microphone openings, are never incontact with liquid as this may damage electronic components.
- While wiping the device with disinfectant, do not press any buttons as long as the device is still wet.
- Wear Cochlear™ transmitter mic within 6-8" of your mouth using either the clip or lanyard.
- Turn Mini Mic 2+ on.
- Clean hands with hand sanitizer after the sanitization process is complete.
- At the end of the class, repeat the sanitization process. Return the Mini Mic 2+ to the student.

